

Argatroban	PROTOCOL FOR ISOLATED HIT, HIT & HITT IN PATIENTS WITH HEPATIC DYSFUNCTION (HEPATIC DISEASE SCORE GREATER THAN 6) NON CATH LAB
Indication	Anticoagulant for prophylaxis or treatment of thrombosis in patients with Isolated HIT, HIT, HITTs (Heparin induced Thrombocytopenia/thrombosis syndrome)
Contraindications	Overt major bleeding; argatroban should be avoided if LFT are greater than or equal to 3 times the upper limit of normal
Baseline Labs	aPTT, PT/INR, Serum Creatinine, CMP, HCT, and platelets
Monitoring	Platelets daily while on argatroban
Monitoring	Guaiac all stools, gastrocult all emesis, visually check for hematuria or other bleeding. CALL PHYSICIAN FOR BLEEDING
Notify Physician immediately for any of the following:	unexplained drop in blood pressure, development of hematoma, drop in hematocrit, significant bleed, flank pain, bright red urine, or bruising
Discontinue	All Heparin, Lovenox (enoxaparin), dalteparin (Fragmin)
Discontinue	Discontinue warfarin until platelets are back to baseline or > 100,000 (preferably when count is 150,000)
Discontinue	IM injections while on Argatroban, obtain physician order for alternative route
Antidote	No antagonist is available to reverse the action of argatroban
Transition to Warfarin	Start when platelet count > 100,000 (preferably when count is 150,000). No starting doses > 5 mg. Overlap with argatroban for at least 5 days.
INR Monitoring	INR daily. If INR > 4 on combination therapy stop argatroban, repeat INR in 4-6 hours, if INR is less than therapeutic restart argatroban at previous rate. Repeat the procedure daily. Discontinue argatroban when two consecutive daily INRs are therapeutic on warfarin alone after therapy has overlapped for at least 5 days.
aPTT Monitoring	Draw aPTT 2 hours and 6 hours after initiation of therapy and every 4 hours after any change in infusion rate until two consecutive aPTT are therapeutic then draw aPTT every 24 Hours
	If aPTT is greater than 100 seconds, stop argatroban, drawn aPTT every 2 hours until aPTT within therapeutic range, restart argatroban at new rate
aPTT Goal	All aPTTs should be run stat
Mixing instructions	aPTT 40-83 seconds (round actual aPTT to closest whole number) Dosage adjustments based on aPTT see Argatroban Rates Change Based on aPTT Chart
Drug Amount	250 mg
Volume	500 ml
Final Concentration	500 mcg/ml
Starting Dose	0.5 mcg/kg/min (Hepatic Disease score greater than 6) Hepatic Disease Score and creatinine clearance calculators are available on the pharmacy web site under calculators If Liver Function Test are greater than 3 times upper limit of normal do not use argatroban
Dosage Range	0.5-10 mcg/kg/min
Maximum Dose	10 mcg/kg/min

Patient Weight (kg)	Patient Weight (lbs)	INITIAL DOSE (mcg/kg/min)														
		0.1	0.2	0.3	0.4	0.5	1	2	3	4	5	6	7	8		
40	88	0.5	1	1.4	1.9	2.4	4.8	10	14	19	24	29	34	38	43	48
45	99	0.5	1.1	1.6	2.2	2.7	5	11	16	22	27	32	38	43	49	54
50	110	0.6	1.2	1.8	2.4	3	6	12	18	24	30	36	42	48	54	60
55	121	0.7	1.3	2	2.6	3.3	7	13	20	26	33	40	46	53	59	66
60	132	0.7	1.4	2.2	2.9	3.6	7	14	22	29	36	43	50	58	65	72
65	143	0.8	1.6	2.3	3.1	3.9	8	16	23	31	39	47	55	62	70	78
70	154	0.8	1.7	2.5	3.4	4.2	8	17	25	34	42	50	59	67	76	84
75	165	0.9	1.8	2.7	3.6	4.5	9	18	27	36	45	54	63	72	81	90
80	176	1	1.9	2.9	3.8	4.8	10	19	29	38	48	58	67	77	86	96
85	187	1	2	3.1	4.1	5	10	20	31	41	51	61	71	82	92	102
90	198	1.1	2.2	3.2	4.3	5	11	22	32	43	54	65	76	86	97	108
95	209	1.1	2.3	3.4	4.6	6	11	23	34	46	57	68	80	91	103	114
100	220	1.2	2.4	3.6	4.8	6	12	24	36	48	60	72	84	96	108	120
105	231	1.3	2.5	3.8	5	6	13	25	38	50	63	76	88	101	113	126
110	242	1.3	2.6	4	5	7	13	26	40	53	66	79	92	106	119	132
115	253	1.4	2.8	4.1	6	7	14	28	41	55	69	83	97	110	124	138
120	264	1.4	2.9	4.3	6	7	14	29	43	58	72	86	101	115	130	144
125	275	1.5	3	4.5	6	8	15	30	45	60	75	90	105	120	135	150
130	286	1.6	3.1	4.7	6	8	16	31	47	62	78	94	109	125	140	156
135	297	1.6	3.2	4.9	6	8	16	32	49	65	81	97	113	130	146	162
140	308	1.7	3.4	5	7	8	17	34	50	67	84	101	118	134	151	168
145	319	1.7	3.5	5	7	9	17	35	52	70	87	104	122	139	157	174
150	330	1.8	3.6	5	7	9	18	36	54	72	90	108	126	144	162	180
155	341	1.9	3.7	6	7	9	19	37	56	74	93	112	130	149	167	186
160	352	1.9	3.8	6	8	10	19	38	58	77	96	115	134	154	173	192
165	363	2	4	6	8	10	20	40	59	79	99	119	139	158	178	198
170	374	2	4.1	6	8	10	20	41	61	82	102	122	143	163	184	204

RATE (ml/hr)

Infusion Rate Change Based on aPTT Please see Rate Change Chart for rate in ml/hr				
aPTT (Seconds)				
33 or less	34-39	40-83	84-96	97 or higher
50% Increase	25% Increase	No Change	25% Decrease	50% Decrease
Rate Change				
If aPTT is greater than 100 seconds, stop argatroban, drawn aPTT every 2 hours until aPTT is within therapeutic range, then restart argatroban at new rate				
Round the aPTT to the closest whole number				